# Ferritin Control Kit

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# Ferritin control kit

**English** For human medical use only!  

for quality control procedures of the quantitative in vitro determination with the Ferritin test kit on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0170	Ferritin test kit	32 tests
Order number: ST0171	Ferritin test kit	16 tests
Order number: ST1700	Ferritin control kit	2 x 1 ml (decision level)

 **Test kit preparation: Allow single test at least 15 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.**

## Summary

This Ferritin control kit is intended for use in the clinical laboratory as a control for quantitative determination of Ferritin. It is prepared from human serum and was found to be nonreactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However, this material should be handled as though capable of transmitting infectious disease.

## Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date on the label. Once opened vials are stable for at least 8 weeks when stored tightly capped at 2 - 8 °C.

## Waste Management

Please refer to local legal requirements.

## Preparation

The control is liquid and ready to use.

## Warnings and Precautions

This test kit is for in vitro diagnostic use only. DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds.

-  Take the necessary precautions for the use of laboratory reagents.
-  Set the sex to „Control“ before running the QC.
-  Compare the results with the values indicated in the LOT-specific value table.
-  If the result is not within the acceptable range repeat measurement.

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**Allow single test at least 15 minutes to warm up to room temperature before use!**

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### ATTENTION!

Allow single test at least 15 minutes to warm up to room temperature before use!

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1. **Preparation of test system**
- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press „Measurement“ button, enter required information using the touch screen

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3. **Processing of control measurement**
- 3.1 Dispense some drops control liquid into serum cup. Aspirate 50 µl from serum cup.
- 3.2 Dispense 50 µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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