HbA1c Control Kit

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HbA1c control kit

English For human medical use only!  

for quality control procedures of the quantitative in vitro determination with the HbA1c test kit on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0110	HbA1c test kit	24 tests
Order number: ST1100	HbA1c control kit	1 x 0.25 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The HbA1c control kit is intended for quality control procedures in examining the precision of the quantitative in vitro determination of haemoglobin A1c with HbA1c test kit on the smart or CUBE laboratory photometer. The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. This product is a control based on human blood material (erythrocytes).

Stability and Storage

This product will be stable until the expiration date when stored unopened at 2 to 8 °C. Once the control is opened it will be stable for 9 months when stored tightly capped at 2 to 8 °C.

Preparation

The HbA1c liquid control is ready to use. Before sampling, gently swirl the vial several times to ensure homogeneity. Observe for inhomogeneity. If necessary gently swirl again.

Procedure

The control should be assayed in the same manner as blood specimens INCLUDING THE LYSATION PROCEDURE. Set sex to „Control“. Use only disposable pipette tips to avoid contamination.

Waste Management

Please refer to local legal requirements.

Warnings and Precautions

This control kit is for in vitro diagnostic use only.

Each individual blood donation used for production was found to be nonreactive when tested with approved methods for HBsAg, anti-HIV 1+2 and anti-HCV. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the control with the same precautions used for patient specimens.

DO NOT INGEST. Avoid contact with skin and eyes.

Take the necessary precautions for the use of laboratory reagents.

Compare the results with the values indicated in the lot specific value table.

If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor. Do not perform patient samples if the controls are outside accepted limits.

Check performance of the test according to national requirements of external interlaboratory tests (EQUAS)

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


Processing of a HbA1c control measurement

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ATTENTION!

Allow single test to warm up at room temperature for a minimum of 10 minutes before use!

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1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place lysis cup, ERS cuvette and ERS cap into test kit rack
- 1.3 Press measurement button, enter required information using the touch screen

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2. Control preparation

- 2.1 Sex-Menu: Select „Control“
- 2.2 Shake vial with control liquid gently 6 to 8 times
- 2.3 Aspirate 10 µl control liquid using the pipette
- 2.4 Dispense 10 µl control liquid into lysis cup (prefilled with lysing buffer)

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3. Control processing

- 3.1 Tightly close and shake lysis cup vigorously 6 to 8 times (= lysate). The lysate may be used instantly after shaking
- 3.2 Pipette 10 µl lysate into ERS cuvette INTO LIQUID
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer. Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: cf967e8ae7c144949ad6b4ec6db8ab63 ; latest revision: DP 7.0 ; date of print: 2019-04-29