PT (INR) Control Kit

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PT (INR) control kit

English For human medical use only!  

for quality control procedures of the quantitative in vitro determination with the PT (INR) test kit on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0180	PT (INR) test kit	32 tests
Order number: ST1800	PT (INR) control kit	2 x 1 ml (N & AK)

Test kit preparation: Allow single test at least 30 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

For use as accuracy control of the PT (INR) test kit. Precision control, i.e. the recording of the day to day variation, is an important procedure in coagulation testing; checking the accuracy of the test system is particularly important. PT (INR) control kit is prepared from selected Austrian citrated plasmas („Austrian Source“) of at least 100 donations. The controls contain stabilizers but no bactericide additives.

Material Required

(not supplied with the control kit)

- Pipettes
- Distilled water

Storage and Handling

The expiry date printed on the labels applies to storage of the unopened bottles at +2 to 8°C.

Stability after reconstitution:

Room temperature:	4 hours
Stored at +2 to 8 °C:	8 hours
Stored at - 20 °C:	8 weeks

The vials can only be frozen once. Upon storage, caps should be screwed tightly.

Waste Management

Please refer to local legal requirements.

Test Preparation and Procedure

Open the vial carefully and reconstitute the contents in 1 ml of distilled water by carefully rotating the vial (avoid frothing!). Allow reconstituted controls to stand for 10 minutes at room temperature before use. Treat the reconstituted control plasma like a citrated patient sample according to the instruction.

Warnings and Precautions

For in vitro diagnostic use only. All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.

Each single donor plasma and each lot of PT (INR) control are tested and found negative for HbSag, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

Set samplotype to „Ctrl.Plasma“ and sex to „Control“ before running the QC sample.

Compare the results with the values indicated in the lot specific value table.

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range please contact your distributor.

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ATTENTION!

Do not wipe the capillary or pipetting tip after aspirating the sample! Allow single test min. 30 minutes to warm up to room temperature (20 - 25 °C) before use!

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1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen display

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3. Test preparation measurement with pipette

- 3.1 Dispense some drops control liquid into serum cup. Aspirate 20 µl from serum cup.
- 3.2 Attention: Entire reagent liquid has to be at the bottom of the cuvette. Dispense 20µl control liquid into ERS cuvette INTO LIQUID (note reconstitution of test kit)! **Carry out the steps 3.3 & 3.4 immediately!**
- 3.3 Apply ERS cap firmly onto ERS cuvette
- 3.4 Place ERS cartridge into analyser. Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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