iFOB Control Kit

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iFOB control kit

English For human medical use only!  

for quality control procedures of the quantitative-immunological in vitro determination with the iFOB test kit on smart or CUBE laboratory photometer.

Order information	Indication	Kit size
Order number: ST0200	iFOB test kit	32 tests
Order number: ST0205	iFOB sample collector kit	32 collectors
Order number: ST2000	iFOB control kit	2 x 2 ml (low/high)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

The iFOB control kit is a solution of human haemoglobin suitable for the accuracy and precision control of the haemoglobin determination in fecal occult blood by turbidimetry methods in clinical laboratories on the smart photometer and CUBE analyser.

Storage and Stability

Unopened vials are stable at 2 - 8 °C until the expiration date on the label. After use please cap vial immediately and put back into refrigerator. Once opened, vials are stable for 30 days when stored tightly capped at 2 - 8 °C.

Waste Management

Please refer to your local legal requirements.

Preparation and Handling

The controls are liquid and ready to use. Thoroughly mix the contents of the vial before each use by gently inverting for several times. After extracting the necessary volume recap and store the vial at 2 - 8 °C immediately.

Control values

Accuracy should be checked according to your local/national QC requirements. The control solution is treated in the same way as patient samples and results should be checked against the acceptable confidence limits supplied on the kit insert sheet and the vial label.

Warnings and Precautions

Set the sex to „Control“ before running the QC.

Compare the results with the values indicated in the lot specific value table.

If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

Contains sodium azide, which may react with lead or copper plumbing to form explosive compounds. Flush drains with copious amounts of water when disposing of this reagent.

The solution was found to be non reactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However this material should be handled as though capable of transmitting infectious disease.

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ATTENTION!

Allow single test to warm up at room temperature for a minimum of 10 minutes before use!

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1. Preparation of test system

- 1.1 Place RFID card
- 1.2 Place ERS cuvette in the test kit rack
- 1.3 Place ERS cap in the test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

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
3. Processing iFOB control measurement with pipette

- 3.1 Aspirate 40µl control liquid from control bottle
- 3.2 Dispense 40µl sample INTO LIQUID in ERS cuvette
- 3.3 Apply ERS Cap firmly on ERS cuvette.
- 3.4 Place ERS cartridge into laboratory photometer. Start automatic sample processing by pressing the  start button on the smart laboratory photometer, or by closing the door of the CUBE laboratory photometer.

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